


Introduction

Tumor and stem cells grow and interact with a tissue microenvironment, which includes the surrounding tissue cells, blood vessels, immune cells, and signaling molecules. The interactions between the tumor/stem cells and the microenvironment play an important role in the regulation of cell growth and differentiation.

Objective & Impact of Research

- The mission of Professor Keyue Shen's Laboratory for Integrative Biosystems Engineering is to develop *in vitro* models of tumor and stem cell microenvironments.
- These models allow researchers to isolate the individual components of the tissue microenvironment and observe their effects on tumor/stem cells. They are developed using microfabrication, which is the process of manufacturing devices that manipulate substances at the μ scale or smaller.
- By utilizing the models, the lab hopes to discover new targets for drug development and develop more effective tumor/stem cell treatments.




Model of microfabricated device designed to replicate hypoxia (lack of oxygen) in tumors. Yellow rectangles represent MCF-7 breast cancer cells. Cells near the edge of the device receive more oxygen than cells near the center.
PC: Yuta Ando [1]

Methods & Results

Analysis of Immunofluorescence Images

- Utilized ImageJ and MATLAB to analyze immunofluorescence images (**Figure 2**) of MCF-7 breast cancer cells stained with hypoxic marker Glut-1 that were grown in mentor's hypoxia microdevice platform (**Figure 1**).
- My program counted 10,018 more objects above a threshold of 165 in the normoxia image vs the hypoxia image (**Figure 3**). Less expression of Glut-1 in normoxia image.
- Mean is slightly higher in hypoxia image due to increased fluorescence in hypoxia image vs normoxia image (**Figure 4**).




PC: Yuta Ando Glut-1

PC: William Koh


PC: William Koh

Food Dye Experiment with Microfluidics Devices


- Conducted experiment with microfluidics devices designed to deliver drugs to cells using food dye (**Figure 5**).
- Control Device: dyes remain separate. Herringbone Device: dyes mix.
- Assessed functionality of the devices by utilizing the Color Histogram tool in ImageJ software.
- Minimal difference in percentage of RGB values in control device shows food dye did not mix (**Figure 6**).
- 19% greater percentage of green pixels in end of herringbone device than start shows food dye did mix (**Figure 7**).



PC: William Koh




PC: William Koh



PC: William Koh

Measuring the Stiffness of Gelatin

- Hydrogels like gelatin are used to mimic biological tissue in microfabricated devices.
- Young's Modulus is a numerical value that describes the stiffness of a material.
- Young's Modulus Formula:
 $E = (F/A)/(\Delta L/L)$
E: Young's Modulus F: Force Applied
A: Area of Section L: Original Length
 ΔL : Change in Length After Force Applied
- Made gelatin in cube mold, placed weight on top of gelatin and calculated E.
- E increases as concentration increases. Higher E indicates greater stiffness.



PC: William Koh

References

[1] Ando, Y., Ta, H. P., Yen, D. P., Lee, S. S., Raola, S., & Shen, K. (2017). A microdevice platform recapitulating hypoxic tumor microenvironments. *Scientific reports*, 7(1), 15233.

Acknowledgements

Special thanks to Professor Shen for this opportunity to experience research. Thank you to my mentor, Dr. Yuta Ando, my lab partner, Yubin Lee, Dr. Katie Mills, and the rest of the SHINE team for supporting me throughout my projects.